TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate www.vevor.com/support

ELECTIRC WHEEL CHAIRS OPERATING INSTRUCTIONS

MODEL: W5216


We continue to be committed to provide you tools with competitive price.

"Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

VEVOR[®]

TOUGH TOOLS, HALF PRICE

MODEL: W5216

NEED HELP? CONTACT US!

Have product questions? Need technical support? Please feel free to contact us:

Technical Support and E-Warranty Certificate
www.vevor.com/support

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

Read and follow all instructions, warnings, and notes in this manual before attempting to operate your power wheelchair for the first time. If there is any information in this manual which you do not understand, or if you require additional assistance for assembly or operation, please contact with your local dealer authorized.

These symbols below in this manual are used to identify warnings and important information. All of them are very important to your safety. It is strongly recommended that you should read and understand them completely.

WARNING! Failure to heed the warnings in the manual may result in personal injury.

ATTENTION! Failure to heed the cautions in the manual may result in damage to the power wheelchair.

1. Introduction

1.1 This product integrates advanced technology with modern style. Especially, its structure for quick dis-assembly and assembly are very convenient for you to store or place at the trunk of your vehicle while traveling. We are certain that the design features, excellent performance and trouble -free operation of this product will ensure your daily life more convenient.

1.2 This product is suitable for all the disabled persons, the elderly and the infirm who have difficulty in walking. The maximum user weight is 120kg.

1.3 This product is of an outdoor type. It is of Model W5216. Here, “W” indicates an outdoor power wheelchair, “5” the code of category of products, “216” the serial number of this type.

1.4 This product is suitable for an outdoor use and the flat streets near buildings for a middle distance, but not for the grass, gravel, or slope of more than 9 degrees and motorized roads, and also not for rainy day.

2. Safety use

2.1 Expiration date

This power wheelchair is energy, non-asepsis, reusable product without chemical substances.

The safe and period of validity for wheelchair on the clinical use not only depends on its structural strength of products, but also on the user usage and use of the product environment, user habits, whether in accordance with the periodic maintenance, maintenance and other factors .

Our products are valid from the date of production for 7 years. (production date on inspection certificate).

2.2 Instructions before use

2.2.1 If not completely read and understand the user manual, please do not drive.

2.2.2 The maximum loading capacity is 120kg, please do not overload driving.

And it is single use, please do not carry passengers.

2.2.3 Please do not drive after drunk or overtired.(Figure 1)

2.2.4 Do not drive at night or in the case of a unclear line of sight.

2.2.5 Please self-test before driving and refer to the manual section 5.3 "Commissioning" and 6.2.1 "Practice before operation " if driving for the first time.

2.2.6 This wheelchair is not waterproof, do not expose in rain and snow as well as driving in rain or snow.(Figure 2)

2.3 Attention during using

2.3.1 This power wheelchair can drive on a good flat surface, do not muddy, rugged, soft, narrow, icy road, driving on bad roads without guardrails dangerous road or waterway. (Figure 3)

2.3.2 This power wheelchair has the certain ability of groove obstacle , at this time you should reduce the speed and travel slowly.

2.3.3 This power wheelchair has the certain climbing ability, the maximum slope is 6°. (Figure 4)

2.3.4 Please avoid driving in the crowds, traffic and other places with heavy traffic.

2.3.5 Please ensure that the controller system fix well, joystick is straight and correct. Sit well and promptly fasten protective tape, after the head as much as possible by the back against the seat back to prevent the bumps on the road, which may could lose control.

2.3.6 Press the on/off button, firstly check if the current fastest speed setting is appropriate for yourself operating proficiency, or it should be adjusted for safety;

We recommend that we should drive slowly at the beginning for every use, and gradually accelerating as well.

ATTENTION!

Please pull the lever up when the power wheelchair stop at the incline, or it will taxiing at low speed, result in personal injury.

WARNING!

If your wheelchair moves accidentally, please immediately release the joystick so that the wheelchair will stop moving automatically. If the joystick is out of control and the brake failer, please cut off the power.

2.4 Electromagnetic interference (EMI)

Your wheelchair may be traveling in the areas affected by electromagnetic interference from some radio transmitters such as radio, wireless intercom, mobile phones and radars etc. In these cases, your wheelchair's driving may be affected by them. Electromagnetic interference may affect the control system of the electric wheelchair. Some disturbances may cause brake failure, automatic start or runaway, and may cause permanent damage to the control system.

2.4.1 Electromagnetic interference classification

- Short-band radio transceiver. These radios have antennas, such as radio waves in the urban bands, walkie-talkies, cell phone network systems and signal transmission devices.
- Medium-band mobile radio transceiver. These are usually installed in the building or outside the vehicle antenna. For example: police, fire, tax, medical ambulance and other radio transceivers.
- Large band radio transceiver. These are usually installed on the tower in the antenna, such as: commercial radio TV transmission system.

ATTENTION!

Wheelchair and battery chargers meet GB / T 18029.21

Wheelchair - Part 21: Requirements and test methods for electromagnetic compatibility of power chair, electric scooters and battery chargers.

ATTENTION!

The chair during driving may be affected by Electro-Magnetic Interference(EMI). Such interference may be from mobile phones, EAS system, etc.

2.4.2 Prevention of Electromagnetic Interference

- Do not use hand-held radio transceivers when the wheelchair power is turned on, such as mobile phones, radios, etc;
- Avoid close to the radio transmission system, such as radio stations, television stations;
- If the electric wheelchair can not control or brake failure, please contact our company or dealer.

3. Structure and performance

3.1The structure of the power wheelchair

This power wheelchair consists of the main part including the front wheel, drive wheel, frame, Joystick, motor with mechanism, armrest, push handle, backrest, seat base, footrest, battery case, the structure shown in Figure 1.



Fig. 1

3.2 Product characteristics of the power wheelchair

- Classification by anti-shock type: wheelchair is internal power supply;
- Classification according to the degree of anti-shock: wheelchair belongs to the B-type application part;
- Classification according to the degree of protection of the feed: IPXO;
- Safety degree when used in flammable anesthesia mixed with air or flammable narcotic with oxygen or oxidized nitrogen Category: Non-AP or APG equipment;
- By operating mode: continuous operation;

- Wheelchair rated voltage and frequency: DC 24V;
- Input power of wheelchair: internal power supply equipment;
- Whether the wheelchair has a protective effect on defibrillation effect: No;
- Whether the wheelchair has a signal output and input part: No;
- Permanent installation of equipment or non-permanent installation of equipment:
mobile devices.

3.3Electrical diagram of power wheelchair

3.4 Product software REV.

Controller version of the power wheelchair W5216 is V1.6.

4. Specification

The main performance parameter of the power wheelchair (size and weight included) Sheet 1.

Max loading weight	120kg
Maximum Speed	≤6km/h
Battery charger	AC 100~240V, 50Hz, 1.5A
Battery	DC 24V 12AH
Front wheel diameter	10 inch
Rear wheel diameter	16 inch
Motor (*2pcs)	DC 24V 250W
Max driving distance	20km
Max safe slope	9°
Static stability	9°
Climbing angle	6°
Seat size	17.5inch

NOTE:

The width of the seat refers to the width between the handlebars.

5.Installation, Folding and Debugging

5.1 Installation



5.1.1 Take out the wheelchair from the packing box and put it on the ground.

Push the chair-frame apart to left and right, and press two tubes on both sides of the seat base down into the grooves on both sides of the chair-frame. (Figure 2)

5.1.2 Unfold battery support(Figure 3), take out the batteries as Figure 4 and then hang them to the supports.



5.1.3 Unscrew the bolts and washers from the controller stem, then put on the controller with the plastic washers inserted, tighten the bolts with screwdriver. (Figure 5)

5.1.4 Connect the plug of the battery box on the frame socket.(Figure 6)

5.1.5 Pull up the handlebars, and lock the folding device.(Figure 7)

5.1.6 Install the footrest on the 2 cylinders of chair-frame(Figure 8), rotate the pedal down(Figure 9), and then the pedals can be adjusted up and down.

5.1.7 Adjust the height of footrest: unscrew the safe nuts from both sides of the footrest pole, take out the bolt(Figure 11), slide the extension tube up and down to adjust the footrest height suitable to the position needed by the user, and then insert the bolt and fix the nuts.

5.2 Folding and storage

Please refer to the below method when transporting or storing the wheelchair:

- ①Put the wheelchair on flat ground and turn off the power supply.

- ②Lay aside the footrests and leg cushion. Turn them upward by 90°.(Figure 9)
- ③Switch off the power supply, take off the batteries(Figure 11), and then put away the battery supports.(Figure 17)
- ④Push the leg pole holder outside, and rotate the leg pole by 90°, then pull upward and take it out.(Figure 18)
- ⑤Adjust the handlebars, turn the lever under the handle to fold the hand shank.(Figure 19)
- ⑥Pull up the seat cushion (Figure 20), then the wheelchair frame could be folded from both sides.

- ⑦Press the handle button for the opening. (Below A)



6. Using and operation

6.1 Using

6.1.1 Controller: an absolutely necessary electrical device. All electronic components to operate the wheelchair are integrated to the device.(Figure 12) Usually, the controller is installed on an armrest, connected together with the motors and batteries to the power box.

- Switch button: Switch button supplies the power for control system for electronic equipment, then equipment supplies the power for motor of wheelchair. Do not use the switch button to stop the chair unless an emergency situation. Otherwise it may shorten the life of the wheelchair drive components.
- USB: For charging daily used small appliances such as mobile phones, small fans.
- Joystick: The joystick is mainly used to control the wheelchair's movement including its speed and direction(forward, backward and left or right etc). The further you push the joystick from its central position, the faster the wheelchair moves. Whenever you release the joystick, it will automatically go back to the center and the brake will be automatically operative to stop the wheelchair.
- Horn button: The horn will sound if you press this button.
- Speed Up/Down Buttons and Speedometer: After turning on the power, the speedometer shows the current maximum speed setting. This maximum speed setting can be adjusted through the speed up button or speed down button by user.

WARNING!

If your wheelchair moves accidentally, please immediately release the joystick and pull up the manual brake.(Figure 13)

●Charger Socket:It is only used to this wheelchair. Do not use the socket to supply power to any other electrical equipment. Otherwise, it may damage the wheelchair's control system or its E.M.C performance (Electromagnetic tolerance).

CAUTION!

If you connect any devices other than the supplied battery charger into the socket, causing damage to the control system, then the warranty is invalid.

6.1.2 Batteries and Its Usage

Fully charge your new battery prior to its initial use. This brings the battery up

to about 90% of its peak performance level; Give the battery another full charge of 8-14 hours and operate your wheelchair again, the batteries will now perform at over 90% of their potential; After four or five charging cycles,

the batteries will top off at 100% charge and last for an extended period.

Please replace a new battery when the battery failure, and the old battery must be returned to supplier because of environmental pollution.

Battery specification	
Type	Deep-cycle lead acid batteries or gel-cell batteries
Dimension	R t12200 (181mm*76mm*170mm)
Power	12V/pc
Amperehour	12Ah

CAUTION!

Do not use your wheelchair when the temperature is less than -25°C or above 50 °C. (Figure 6)

6.1.3 Charging batteries

The battery charger is an important part of the wheelchair. The off-board charger attached to this product can charge its batteries quickly and easily to make your wheelchair the best.

CAUTION!

The charge of 24V/2A supplied by our company meets the requirements mentioned in Section I (《General Safety Requirements》) of medical electrical equipment GB9706.1- 2007. The

battery must be charged by the charger supplied by us. Any other types of chargers are forbidden.

Please note below tips when charging:

- Be certain the controller is powered off and the wheelchair is in the drive mode, instead of the manual mode.

Connect the 3-pin output plug of the charger to the controller. (Figure 15)
The green LED on the charger lights when the batteries are fully charged.
We recommend you that the batteries is charged for 8-12 hours.

- Remove the charger and power plug when fully charged, and put them into the bag behind the seat rest.

WARNING!

Never get the batteries frozen, neither charge the frozen batteries, otherwise it may lead to personal injury and damage to batteries.

6.1.4 Overload protector

The overload protector is set for the wheelchair's safety. Once motors overload, the overload protector will trip to cut off the power to prevent the motors and electric components from damage. In order to recover the function of protector, professionals only to unscrew the bolts from the battery box, replace the backup fuse with the broken one, then cover the battery box and fasten the bolts. (Figure 16)


6.1.5 Brake lever

Pull back the brake lever when no need to drive the wheelchair. Press the wheel with the brake lever to keep the wheelchair remain still.

Push the brake lever forward to keep it away from wheel when need to drive the wheelchair.

6.1.6 Safety belt

For your safety, the safety belt must be fastened before you drive the wheelchair.(Figure 14)

6.2 Operation

6.2.1 Preparation before operation

- Do not turn on the power switch when sit on the power wheelchair, and do turn off the power switch when get off from the wheelchair. Please observe whether the hand brake works or not.
- Please fold the pedal first then hold the armrest to sit on the wheelchair; Do not get off the wheelchair by stepping on the pedal. Otherwise the chair maybe turned over, which is dangerous.


6.2.2 Practice before operation

- Find a spacious place like square and have a assistant to help you practice until you have enough confidence to operate it.
- Be certain to shut down the power when you up and down the wheelchair, and set a speed you need through speed adjustment button.
- We recommend you set the lowest speed until you can operate the electric wheelchair skillfully. (Figure 7)

- To Practice the stop operation, forward and backward. Push the handle, the wheelchair move to anywhere you want.(Figure 8)
- Firstly, ensure that it is the lowest speed when practice forward operation. After Skilled, you can practice “S” shaped turn. After you are familiar with "S" shaped turn, to practice backward operation, and pay attention to the speed control setting. And the backward speed should lower than forward speed.(Figure 9)

fig. 7

start at the
lowest speed



fig. 8

fig. 9

6.2.3 Operation

Please do it according to precautions operation.

7.Fault diagnosis and troubleshooting

This product has an automatic fault warning function to give you very

convenience. Once it malfunctions, the LEDs on the controller will flash with alarm sounded by the horn. You can find where a fault is according to information in Table 2.

If the fault still exists after checking out the fault based on that in table 2, Please consult your service agent immediately.

2 alarm+controller panel pattern 1	Check the communication connection cable between upper controller and lower controller
2 alarm+controller panel pattern 2	Check the left motor and wiring harness
2 alarm+controller panel pattern 3	Check the left brake and wiring harness
2 alarm+controller panel pattern 4	Check the right motor and wiring harness
2 alarm+controller panel pattern 5	Check the right motor and wiring harness
2 alarm+controller panel pattern 6	The controller is in the over current state or the motor end short circuit , the load is too heavy
2 alarm+controller panel pattern 7	The joystick lever of the controller is malfunctioning or it is no longer in position/ in contact
2 alarm+controller panel pattern 8	Failure of the controller itself
2 alarm+controller panel pattern 9	Battery under voltage or battery failure

8.Safety devices and accident treatment

8.1 Overload Protector: Refer to instruction manual 6.1.4

8.2 Safety belt: Refer to instruction manual 6.1.6

8.3 Please release the joystick whenever there is a wheelchair fault , the wheelchair will be stopped.

9.Maintenance and repair

9.1 Maintenance

Like other motorized vehicle, your power wheelchair also requires routine maintenance. Some checks can be performed by yourself, others you can ask for assistance from your service agent. Preventive maintenance is very important. If you follow the maintenance and checks in this section, your wheelchair will give you years of trouble-free operation. If you have any doubt your wheelchair's care or operation, please contact your service agent or our after-sale service dept.

9.1.1 Humidity

Your wheelchair, like most electrical and mechanical equipment, is susceptible to external conditions. In any case, the wheelchair should be avoided damp environment. Direct or prolonged exposure to water or dampness could cause the wheelchair to malfunction electronically and mechanically. Water can cause electrical components and the chair's frame to corrode.

9.1.2 Temperature

- Some parts of your wheelchair are susceptible to temperature.
- In extremely cold temperature, the battery may be frozen. Special temperatures may cause a lot of factors to freeze, like charger type, usage, battery components (such as sealed lead-acid batteries or gel batteries);
- Temperature above 55° may cause your wheelchair's speed to reduce.

9.1.3 General Guidelines

- Avoid beating the controller, especially the joystick.
- Avoid prolonged exposure of your wheelchair to extreme conditions, such as hot, cold or moisture environment.
- Keep the controller clean.

- Check all electric connections, including the cable and connectors of the charger, and ensure that they are all tight and secure.
- If only red LEDs on the Battery Gauge lights, the batteries are nearly running out of charge. You should recharge the batteries as soon as possible. We recommend charging the battery for 8-12 hours
- Check the rear wheel inflatable situation, if the tire deformed serious when riding , it should be inflated;
- The frame surface has been sprayed with a clear sealant coating. You can apply a light coat of car wax to make the surface keep a high gloss.
- Check all cable connections. Make ensure they are fastened and are not corroded. The battery must be placed in the battery area, flat, the battery plug inward, relative display, refer to the correct connection on the main frame label;
- All wheel bearings are lubricated and sealed. Do not need to lubricate them.
- Check if there is loose phenomenon for wheel hub, drive device, and chair itself, if loose, please screw tightly in time.

WARNING!

Please be sure the tire pressure is within $230\text{KPa} \pm 10\%$.

Please inflate the tires with regulated air source and check with pressure gauge. Do not exceed the pressure range, otherwise may cause exploded tires or personal injury. (Figure 11)

9.1.4 Maintenance after use

- Turn off the power.(please disconnect all the connectors if not use for long time.)
- Inhibit children and unconscious person to use the wheelchair.
- Store the wheelchair in normal temperature to prevent reformation so that it keeps its performances for long period.
- Clean the wheelchair with a clear and soft cloth and dry it. Never use any chemicals to clean it.(to prevent deformation and discoloration)

- Remove the cloth cover of the seat rest to be washed if it is dirty, and then dry it for use.

9.1.5 Daily Checks

In order to keep the wheelchair in good condition, you should check before each use. And to maintain the weekly, monthly, semi-annual investigation, check the project see (Table III).

9.1.6 Cleaning

Way of cleaning as below:

- Never wash your wheelchair with water or expose directly to water.
- Surface of wheelchair frame is coated with a protective coating. Therefore, it is very easy to wipe it clean with a damp cloth. Never use any chemicals to clean the vinylon seat and armrest, as they may cause the latter slip or chapped. You can use a damp cloth and neutral soapy water to clean them, and then dry thoroughly.

9.2 Repairment

9.2.1 Users could repair or replace some spare parts like push handle, push handle cover which could be self-maintained or replaced by users under the instructions of after-sales service staff of our corporation or the dealers'.

9.2.2 Within the warranty, our corporation or the dealers will maintain or replace the faulted spare parts, like motor, battery, charger, controller, bearing, shaft sleeve, etc., due to material or manufacturing defects free of charge.

9.2.3 The spare parts like motor, controller, etc. which are confirmed by the after-sale staff or our corporation or our dealers' to be disassembled could be sent to our corporation or the dealers for maintenance.

9.2.4 Please contact your distributor or the after-sale service dept. of our corporation for any questions regarding the maintenance of power wheelchairs.

Inspection Items	At any time	Weekly	Monthly	Six monthly
Wheel hubs, Head stock Driving mechanism and fasteners			√	
Joystick function	√			
Brake System	√			
Connection		√		
Battery Condition	√			
Tire Condition			√	
Inflation condition of rear wheels	√			
Frame Condition				√
Motor and actuator condition				√
Front Wheel Condition		√		
Pureness	√			

Clean the wheelchair
with dry cloth.

Caution!

In order to guarantee the quality of the power wheelchair, Please contact us or our dealers to replace the spare parts with same models and specifications if the consumable parts like batteries, tyres, controller, motors and charger, etc. need to be replaced. The maintenance should be proceeded at our corporation or the dealers' maintenance dept.

10. Transportation and storage

10.1 Transportation

You can load and transport according to the shipping marks and graphics. For details, see the attachment.

10.2 Storage

Your wheelchair should be stored in a clean indoor environment with relative humidity $\leq 80\%$, good ventilation and free from corrosive gas. Remove the batteries from the wheelchair prior to storage. Otherwise, the frame may rust and the electronics may be damaged.

11. Opening and checking

The packing list is attached in the package. Please check if any parts missed or damaged.

12. Quality Assurance

12.1 Warranty

12.2.1 Chair frame supply 3 years warranty service.

12.2.2 Within 1 year from the purchase date, for the following parts, we will supply free maintenance and replace service for original customer after the dealer check there are materials and production defect.

- Electrical control or lever system
- Motor/Drive System
- Bearing and shaft sleeve



12.2.3 Six-Months warranty for the batteries.

Please note the warranty service is provided by your dealer, and finally finished by our after-sales department and dealer together.

Out of warranty

- ABS plastic cover shell and rubber pad
- Tyre
- Interior decoration
- Damage by abuse, wrong operation, accidents and negligence
- Damage by wrong operation, maintenance and storage.
- Business use or other abnormal use

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Pooledas Group Ltd
Unit 5 Albert Edward House, The Pavilions
Preston, United Kingdom

Made In China

VEVOR[®]

TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate

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Importerad till AUS: SIHAO PTY LTD . 1 ROKEVA STREETEASTWOOD
Importerad till USA: Sanven Technology Ltd. Suite 250, 9166 Anaheim
SHUNSHUN GmbH