VEVOR[®]

TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate www.vevor.com/support

LAB CHILLER CIRCULATOR USER MANUAL

We continue to be committed to provide you tools with competitive price. "Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

VEVOR[®]

TOUGH TOOLS, HALF PRICE

LAB CHILLER CIRCULATOR

MODEL: DC-0506/DC-2006

Photo for reference

NEED HELP? CONTACT US!

Have product questions? Need technical support? Please feel free to contact us:

Technical Support and E-Warranty Certificate
www.vevor.com/support

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.


SAFETY INSTRUCTION AND PRECAUTIONS

WARNING:

1. Before using, the slot should be in a liquid medium. (note: liquid not contain acid alkali).
2. The total power of the power is to be greater than or equal to the instrument. The power must have a good "grounding" device. Using the voltage will be subject to technical parameter table.
3. The device should be placed in a ventilated dry place; the back and sides leave obstacles within 300 mm distance.
4. After use, switch off and unplug the power plug.

Instrument appearance introduction

Features:

- 1 Refrigeration System: low noise air-cooled totally enclosed compressor unit is

adopted, which has the advantages of fast refrigeration speed and strong temperature stability.

2 Control System: Latest proprietary software research and development of homemade temperature control technology, with high standards of PT100 and all imports of electronic parts, etc.

3 Circulating System: it is composed of internal and external circulation. During external circulation, the constant temperature liquid in the tank can be introduced outside, the second constant temperature field can be established, and it can also be used as a cold source to lead the cooled liquid in the tank to the experimental container outside the machine for cooling; During internal circulation, the liquid temperature in the tank can be very uniform and stable.

4 Protection System: Refrigeration is overheated, over-current protection, call control system has over temperature alarm, can be set over temperature alarm temperature, the temperature of the upper and lower limits can cut off the load automatically when over temperature, and other functions.

5 High-end LCD control instrument, simple operation, temperature stability is strong, P.i.d can automatically search or artificial adjustment

6 The accurate temperature correction function can reach 0.01°C , and the maximum temperature fluctuation can reach $\pm 0.02\text{-}0.05^{\circ}\text{C}$ (depending on the model).

Technique Parameter:

Model	Temperature range ($^{\circ}\text{C}$)	Temperature fluctuation ($\pm^{\circ}\text{C}$)	Slot opening (mm)	Groove depth(mm)	Power (KW)
DC-0506	-5-100	0.03	180x160	150	0.8~1
DC-2006	-20-100	0.03	180*160	150	0.9~1

Operating Steps:

1.Connect the circulating pump before startup (if there is no water inlet and outlet, this step can be omitted):

1.1 for the connection of circulating pump during internal circulation, connect the water outlet and water inlet directly with both ends of the hose supplied with the goods.

1.2 for the connection of circulating pump during external circulation, connect the water outlet to the inlet of container / equipment outside the tank with hose, and connect the water inlet to the outlet of container / equipment outside the tank.

2. Add the corresponding liquid medium into the tank. The liquid level of the liquid medium shall not be 20mm lower than the worktable to avoid damage caused by exposed dry burning of the heater.

2.1 When the operating temperature is below 8°C, industrial alcohol or antifreeze is generally used for liquid media.


2.2 At an operating temperature of 8°C-75°C, pure water is generally selected for liquid media.

2.3 At an operating temperature of 75°C-100°C, dimethyl silicone oil with a viscosity of 5cs is generally selected for liquid media.

3. Plug in the power, turn on the "main power" switch, and then turn on the "power" button on the operation panel. (Note: "main power" switch is on the back of the instrument).

4. Perform the instrument operation as follows:

4.1 Display Description:

4.2 description of instrument keys:

Set function key shift key ▼ decrease key ▲ increase key

4.3 setting of temperature parameters:

Press the setting function key "set" to enter the setting state of the temperature setting value. At this time, the last value of the setting temperature "SV" on the display screen flashes. First press the shift key, then press the plus or minus key to set the working temperature you need, and then press the setting function key "set" once. At this time, the last value of the timed shutdown time (min) "0000" on the upper right corner of the display screen flashes. If it is not necessary to keep the default "0000" regularly, Press the setting function key "set" again to exit the parameter setting and save the parameter setting value(PV). At this time, the display "PV" shows the temperature of the liquid medium in the current tank.

5. After the temperature parameter setting is completed, first open the "circulate" button on the operation panel, and then open the "refrigeration" button (Note: when the working temperature exceeds 40 °C, the compressor automatically enters the protection, and the refrigeration does not start), and then the microcomputer enters the automatic control state.

Packing list

1. Instructions * 1
2. Device host *1
3. Silicone tube * 1
4. Insulation pipe * 2
5. Flow control valve * 1
6. Fuses * 2
7. Cover plate * 1
8. Power cord * 1

CORRECT DISPOSAL


This product is subject to the provision of European Directive 2012/19/EU. The symbol showing a wheeled bin crossed through indicates that the product requires separate refuse collection in the European Union. This applies to the product and all accessories marked with this symbol. Products marked as such may not be discarded with normal domestic waste, but must be taken to a collection point for recycling electrical and electronic devices.

Address: Shuangchenglu 803nong11hao1602A-1609shi, baoshanqu, shanghai 200000 CN.

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YH CONSULTING LIMITED.
C/O YH Consulting Limited Office 147, Centurion House,
London Road, Staines-upon-Thames, Surrey, TW18 4AX

VEVOR[®]
TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate
www.vevor.com/support

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Importerad till AUS: SIHAO PTY LTD. 1 ROKEVA STREET EASTWOOD NSW 2122 Australien
Importerad till USA: Sanven Technology Ltd. Suite 250, 9166 Anaheim Place, Rancho Cucamonga, CA 91730
E-CrossStu GmbH