VEVOR[®]

TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate www.vevor.com/support

LIQUID NITROGEN CONTAINER USER MANUAL

We continue to be committed to provide you tools with competitive price.

"Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

MODEL: YDS-3/YDS-6/YDS-10/YDS-20/YDS-30

Photo for reference

NEED HELP? CONTACT US!

Have product questions? Need technical support? Please feel free to contact us:

Technical Support and E-Warranty Certificate
www.vevor.com/support

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

SAFETY INSTRUCTION AND PRECAUTIONS


WARNING:

1. The liquid nitrogen tank should be placed in a well ventilated and cool place, avoiding direct sunlight
2. It is strictly prohibited to use in rainy or humid environments to avoid frosting and freezing of neck tubes, plugs, and other accessories;
3. It is strictly prohibited to tilt, lay horizontally, invert, pile up, collide and other behaviors. During use, it should be handled gently and always kept upright
4. Do not open the vacuum nozzle of the container. Once the vacuum nozzle is damaged, it will immediately cause vacuum failure;
5. Use professional protective equipment. When operating liquid nitrogen, antifreeze gloves, goggles, and protective clothing should be worn to prevent frostbite.

Product brief instruction:

Used for storing semen/embryos/cells/dry ice/metal materials in laboratory, medical, farm, beauty and other fields. The liquid nitrogen biological containers produced by our company have the characteristics of reasonable structure, easy portability, and good refrigeration performance. They can be used in animal husbandry, medical research departments for long-term low-temperature storage and transportation of livestock, frozen semen, vaccines, cells, microorganisms, etc. using liquid nitrogen (-196 °C) as refrigerant. They can also be applied in industries such as national defense, scientific research, machinery, electronics, metallurgy, energy, catering, beauty, etc.

Instruction for the Container Model

B".Diameter (mm) ,size is not marked it is 50mm.

"C". Nominal volume of container


"D" means nitrogen

"S" means biological

"Y" means liquid

For example: YDS-30 means volume 30L, Dia.50mm, liquid nitrogen biological container.

Product structure

1. Handle
2. Plug
3. neck tube
4. Adsorbent
- 5, outer liner
- 6, multi-layer insulation
- 7, inner liner
- 8, bucket
- 9, vacuum head
- 10, base

1. As shown in the figure, the container consists of an inner liner, an outer liner, a bucket, and a neck tube connecting the inner and outer liners. The inner and outer

liners are made of high-strength aviation aluminum alloy, which is sturdy and durable,

2. There is a high vacuum between the inner and outer tanks, and the vacuum film is equipped with multiple layers of imported insulation materials, thus having good refrigeration performance

3. The neck plug is made of non-toxic polyurethane with good insulation. It has a dual function of reducing liquid nitrogen evaporation and fixing the lifting cylinder.

4. The neck tube is made of glass reinforced compounds, which have extremely low thermal conductivity and can control the external heat entering from the neck tube to a minimum.

5. Gas adsorbents are installed in the vacuum area, which can adsorb gases emitted by metals or other substances in the true cavity, thereby maintaining a high vacuum for a long time.

6. Place the container in liquid nitrogen to store frozen items. The upper part is fixed by the small groove of the neck plug and the notch of the scoring plate, ensuring stable and reliable transportation

7. Bucket needs to be rotated 180 degrees to be drawn out.

Liquid nitrogen Container model

Product model	Capacity L	Caliber mm	Outer diameter mm	Height mm	Static liquid nitrogen storage period	Bucket		
						Qty	Outer diameter	Height
			≤			≥	pcs	mm
YDS-3	3	50	224	450	30	6	36	120
YDS-6	6	50	304	500	60	6	36	120
YDS-10	10	50	304	550	100	6	36	120
YDS-20	20	50	384	640	175	6	38	276
YDS-30	30	50	438	700	270	6	38	276

Note 1: The static liquid nitrogen storage period of the product is tested using the weighing method.

Note 2: In the product model column, B represents a container with support at the bottom or side of the inner liner, which is resistant to tipping and transportation

Note 3: "-" in the delivery column indicates that the product can be delivered according to user requirements. Please confirm and indicate whether to simplify when signing the order contract with the user

Note 4: Due to technological improvements, some data may have some modifications.

Usage notice

Before using the product, users should open the box and inspect the appearance quality of the product. If any defects are found in the appearance of the container, they should return this information to our company's sales department (contact information can be found on the back cover) within 5 working days after receiving the product. The agent and distributor units should follow the relevant terms of the agreement.

1. The container can only be used to fill liquid nitrogen, and it is not recommended to fill other liquid gases or liquid oxygen.
2. Liquid nitrogen is an ultra-low temperature liquid (-196 °C) that can cause skin frostbite if splashed onto the skin. Therefore, when filling and removing liquid nitrogen, warm gloves should be worn and not bare feet or slippers to avoid splashing and injuring people.
3. Only the neck plug that matches the factory can be used to seal the tank mouth, and other plugs cannot be used as substitutes, let alone sealed plugs, to avoid damage to the container caused by the continuous evaporation of liquid nitrogen and the increase in nitrogen pressure.
4. When checking the liquid level inside the container, a solid plastic or wooden rod should be inserted into the bottom and taken out after 5-10 seconds. The

length of frost formation is the liquid level height, and hollow rods or plastic tubes should not be used.

5. Although nitrogen is a odorless, odorless, and non-toxic gas, in poorly ventilated conditions, the evaporated nitrogen can cause an increase in the amount of nitrogen in indoor air and a relative decrease in oxygen content, posing a risk of suffocation to the human body. Therefore, in the storage and operation points of liquid nitrogen biological containers, indoor ventilation should be strengthened to maintain fresh air.

6. Always operate the tank will cause the storage time reduce faster.

7. The neck tube is the channel for the inlet and outlet of liquid nitrogen and frozen items, and should be used with caution. If the neck tube breaks, the container cannot be repaired.

8. The vacuum sealing joint of the container is a key component to maintain the vacuum of the container inter layer. Once damaged, it will immediately cause the vacuum of the container inter layer to deteriorate and the product cannot continue to be used. Therefore, users are not allowed to open the vacuum sealing joint of the container without authorization.

Maintenance and upkeep

1. When putting in and taking out frozen items, try to keep the opening time of the can as short as possible to reduce liquid nitrogen loss, and do not completely lift out the lifting cylinder to avoid affecting the storage effect of the items.

2. Strictly prevent excessive impact and collision. Due to the high vacuum state of the tank, the outer liner will withstand the contrast of the surrounding atmosphere. Therefore, the container (outer liner) will not produce a depression during excessive impact and collision. If the evaporation does not increase significantly after the depression, it can still be used. If the evaporation performance deteriorates, the container cannot be repaired.

3. Regularly check the weight of the container and promptly understand how many days liquid nitrogen can still be used. When the amount of liquid nitrogen is less than one-third of the total capacity, replenish it immediately to prevent the storage items from necrotizing.

4. During the use of the container, due to the impurity of liquid nitrogen, water will gradually accumulate in the inner tank. At the same time, the multiplication of

bacteria will reduce the fertilization rate of semen. Therefore, it is recommended to wash the container with warm water at 40-50°C 1-2 times a year.

Liquid nitrogen static daily evaporation test

1. The static liquid nitrogen storage period of the product is tested using the weighing method.
2. The test should be conducted in an ambient temperature of 20 °C (± 5 °C), and cannot be conducted in a directly ventilated area.
3. The test should fill more than half of the container with liquid nitrogen, let it stand for 48 hours, and weigh (A) grams for the first time. After 3 days, weigh (B) grams for the first time. Calculate the evaporation rate Q according to the following formula:

$$Q=(A-B)/3 \text{ (grams/day)}$$

If it needs to be expressed in terms of volume (liters), use the following formula to convert: $D=Q/808$ (liters/day)

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YH CONSULTING LIMITED.
C/O YH Consulting Limited Office 147, Centurion House,
London Road, Staines-upon-Thames, Surrey, TW18 4AX

Technical Support and E-Warranty Certificate
www.vevor.com/support

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$$Q=(AB)/3 \text{ (grammi/giorno)}$$
Importato in AUS: SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD NSW 2122 Australia
Importato negli Stati Uniti: Sanven Technology Ltd. Suite 250, 9166 Anaheim
E-CrossStu GmbH