# **VEVOR<sup>®</sup>**

**TOUGH TOOLS, HALF PRICE**

Technical Support and E-Warranty Certificate [www.vevor.com/support](http://www.vevor.com/support)

## **MICRO INLAID MIRROR**

**MODEL:ZQ-1**

We continue to be committed to provide you tools with competitive price.

"Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

# VEVOR<sup>®</sup>

TOUGH TOOLS, HALF PRICE

## MICRO INLAID MIRROR

**Note:** The actual product you receive determines its look.

### **NEED HELP? CONTACT US!**

Have product questions? Need technical support? Please feel free to contact us:

**Technical Support and E-Warranty Certificate**  
[www.vevor.com/support](http://www.vevor.com/support)

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

# Important safety instructions


Warning - To reduce the risk of injury, user must read instructions manual carefully.

## **!WARNING!**

1. Use the product only after receiving professional metalworking training.
2. Wear necessary personal protective equipment such as safety goggles, earplugs, and protective gloves during use.
3. Ensure that the workspace remains clean, dry, and well-ventilated during operation.
4. Keep away from children and read this manual in detail before use.

## Save this manual


## 2. Install Ring /Bracket Assembly

- (Fig.1)Insert the rod into the bottom of the mounting post,then attach handle.
- Make sure the washer is under handle.
- Position handle away from stand arm by lifting the handle while pressing the center button

(Fig.2).This does not loosen the handle.To attach the stand to your work bench,see microscope Stand Base Template (page 3).The microscope body sets into the ring and is secured by three thumb screws,included.

### 3. Install Headrest

After installing the microscope,attach the headrest (Fig.3).Make sure white nylon washer is between the T-Knob and the headrest mount.Adjust the angle using the side T-Knob;adjust depth using the top t-knob.The forehead end may also be adjusted

### 4. Setup and Operation

After the stand is secured to the bench top there are many choices for adjustments

to fit your needs. Turn tension knobs clockwise to add friction, counter-clockwise to reduce friction.

A. Controls front arm. Tension knob.

B. Controls friction during movement of back arm. Tension knob.

C. Controls angle of microscope on pivot ball. Loosen to adjust.

D. Adjusts focus. Tension knob.

E. Locks arm height on the stand rod.


F. Controls the drag of the arm swing.

G. Tighten for additional arm rigidity.

## 5. For Wood Top Benches

Decide where to locate the microscope stand on your bench top. Make sure there is plenty of "reach" so the microscope is over your work. Using the template below, tape it exactly where the base will be located. With a center punch or awl mark the centers of the screw holes through the paper and into the bench top. Remove the template and drill a starter hole for the screws using a 1/8" drill bit. Drill all four holes a little more than 1" deep. Hold the drill as square to the surface as possible. Attach the stand to the bench with the 1-1/2" wood screws that are included in the package. Be sure to tighten the screws so that the heads rest tight against the metal base. If you are having trouble setting the screws you probably did not drill the holes deep enough.

## 6. For Metal Top Benches

it is advisable to use 10-32 flat head bolts (not included)with a length of about 1/2"more than the thickness of your bench top.Use the template below as described above to locate the centers and drill 1/4"holes through the bench top.Bolt the stand down using washers and nuts on the underside.If the bench top is made of a sheet metal with a lightweight substrate we recommend adding a board or metal plate to the underside for extra strength.

## 7. Adjust the specimen slide

- 1)Turn the zoom control knob to the maximum magnification.
- 2)Turn the diopter adjusting rings to the zero.
- 3)Observe the specimen through the right eyepiece and make the image clear by turning the focusing knob.

Fig.4 5)Rotate the zoom control knob to the minimum magnification.

5)Observe the specimen through the right eyepiece and make the image clear by turning the right diopter adjusting ring②.(Fig.3)


6)Redo the step(1),(3),(4)and (5)till the right adjusting ring is more precise.

7)Do the step (4)and make the image clear which is observed through the left eyepiece by turning the left diopter adjusting ring①.(Fig.4)


Note:User could adjust the height distance between the micro-scope and the object to about 165mm to get a better clear picture.(If it's not clear,adjust the height according to your multiplier need.)

## 8. Adjust the interpupillary distance

- 1)Adjust the prism housing along the direction of arrow-head of the Fig.5 till the observation is comfortable.

**Fig. 4**

**Fig.5**

## 9. Parameter

1. Name: Micro inlaid mirror
2. Model: ZQ-1
3. Continuous zoom: 7X-45X
4. Machine objective lens: 0.5X/165mm
5. Machine eyepiece: WF10X/20mm

**Address:** Shuangchenglu 803nong11hao1602A-1609shi, baoshanqu, shanghai 200000 CN.

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YH CONSULTING LIMITED.  
C/O YH Consulting Limited Office 147, Centurion House,  
London Road. Staines-upon-Thames. Surrey. TW18 4AX

**VEVOR<sup>®</sup>**

**TOUGH TOOLS, HALF PRICE**

**Technical Support and E-Warranty Certificate**

**[www.vevor.com/support](http://www.vevor.com/support)**

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Importado a Australia: SIHAO PTY LTD, 1 ROKEVA STREETEASTWOOD
Importado a EE. UU.: Sanven Technology Ltd., Suite 250, 9166 Anaheim
E-CrossStu GmbH
C/O YH Consulting Limited Oficina 147, Centurion House,