VEVOR[®]

TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate

www.vevor.com/support

VISCOMETER

MODEL:NDJ-5S/NDJ-9S

We continue to be committed to provide you tools with competitive price.

"Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

VEVOR[®]

TOUGH TOOLS, HALF PRICE

VISCOMETER

Model:NDJ-5S/NDJ-9S

NEED HELP? CONTACT US!

Have product questions? Need technical support? Please feel free to contact us:

Technical Support and E-Warranty Certificate
www.vevor.com/support

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

SAFETY INSTRUCTIONS

1. Before application, it is necessary to understand all safety procedures including instrument power, water and pressure.
2. Ensure that all safe operating procedures are fully understood and followed: including personal safety, but also the safety of instruments and samples.
3. No strong electromagnetic interference, drastic vibrations, corrosive gases in the vicinity of the instrument.

Correct Disposal:


This product is subject to the provision of European Directive 2012/19/EU. The symbol showing a wheelie bin crossed through indicates that the product requires separate refuse collection in the European Union. This applies to the product and all accessories marked with this symbol. Products marked as such may not be discarded with normal domestic waste, but must be taken to collection point for recycling electrical and electronic devices.

MODEL AND PARAMETERS

Model	NDJ-5S	NDJ-9S
Adapter Input	AC 100-240V 50/60Hz	
Adapter Power	9W	
Rating(adapter	DC12V 2A	
Measuring Range	1~100000	10~6000000mpa.s
Measurement	+/-1% (Newtonian liquid)	
Rotor Type	1#, 2#, 3#, 4#(Spindle 0 is optional.)	
Rotational Speed	6/12/30/60RPM	0.1/0.3/0.6/1.5/3/6/12/30/60
Operation	Temperature: 5-35 °C; Humidity: ≤80%	

STRUCTURE DIAGRAM

- 1.Level Vial 2.LCD Screen 3.Operation Panel
4.Outer Cover 5.Spindle Protection 6.Spindle
7.Spindle Connector 8.Base 9.Base Adjusting Knob

INSTALLATION

- 1). Take out base, lifting column, handle from the packing box. Screw the column into the base (rack side towards the user). Fix handle onto the lifting block.
- 2). Adjust lifting button for tightness. Put T-shaped block on the rear of the instrument into handle and get it fixed, making the instrument balanced.
- 3). Adjust 3 screws on the base to keep the level bubble at the center of the black circle.
- 4). Take off the protection cap at the bottom of the instrument.
- 5). Plug in.

HOW TO USE

- 1.Prepare liquid sample to be tested, pour it into beaker or flat container with diameter no less than 70mm. Control the temperature of the liquid sample correctly. (Variations of temperature will affect viscosity directly.)

2. Screw anti-clockwise protection bracket (T-shape) into the bottom end of the instrument.
3. Choose appropriate spindle and screw it anti-clockwise into the connector.
ATTENTION: slightly lift direction connector up when spindle is installed to avoid damages to the pivot.
4. Turn lifting button to slowly immerse spindle into the liquid until the marker on the spindle (groove or marking line) is aligned with the liquid level.
5. Adjust the balance of the instrument again.
6. Should keep testing temperature stable to ensure accurate read out. Otherwise readout may fluctuate.

OPERATION

1. Date and Time Setting:

Turn on, the instrument will display the date and time for 3 seconds, then the instrument automatically enter the state of measurement. During this period, Press OK can be entered into the measurement state. If press the AUTO, you can set date and time.

Such as:
2018-01-01
12:00:00

2. Parameter Setting:

Hold the OK button and turn on the power, the instrument shows:

→Zero: 0.0%
Buzzer: ON
Bit rate: 9600

Language: english

Press ▲ and ▼ to select, then press AUTO. Press ▲ and ▼ to modify parameters. Press OK to confirm.

There are two ways to exit the set state and enter the state of measurement.

A. Hold the AUTO not loose, about 2 seconds, it will exit setting, and the modified parameters are not saved.

B. Hold the OK not loose, about 2 seconds, it will exit setting and keep the modified parameters.

1). Zero setting range: 0.0%-1.0%, for example, setting is 0.5% when measuring, if the torque is less than 0.5%, and the viscosity values show 0.

2). With ▲ and ▼ to choose Zero, press AUTO, “XX%” flicker, then press ▼ and ▼ to modify %, and press OK to confirm.

3). Buzzer setting

With ▲ and ▼ to choose Buzzer. Press AUTO, “ON” And “OFF” flicker, and press OK to confirm.

4). Bit Rate setting

With ▲ and ▼ to choose bit rate. Press AUTO, select 1200, 2400, 4800, 9600, then press OK to confirm.

5). Language

With ▲ and ▼ to choose chinese or english.

3. Correction of Viscosity Deviation

Hold ▲ and ▼ at the same time, and turn on the power, the viscometer shows:

```
viscosity correction
spindle: 1  speed: 6
→CP0:      X X.X%
CPX:       XX .X%
```

The range of viscosity correction percentage is $\pm 12.5\%$. CP0 is the percentage of the idling. CPX is the percentage of viscosity deviation.

Such as, the viscosity is 376.2mpa.s, No. 1 spindle, speed 6RPM, the test

value is 360mpa.s. The error is $(360-376.2) / 360 = -4.5\%$, the CPX should be -4.5%. If the test value is 390mpa.s, the error is $(390-376.2) / 390 = 3.5\%$, the CPX should be 3.5%.

How to revise CPO and CPX:

- 1). Choose the same spindle and speed of testing.
- 2). Press ▼ to choose CPO or CPX.
- 3). Press AUTO, ▲ and ▼ to modify the value, and press OK to confirm.

Exit Correction

- 1). Hold AUTO for 2 second, exit the modified parameters are not preserved.
 - 2). Hold OK for 2 seconds, exit and keep the modified parameters.
- After modify the CPX, under the idlling, if the torque % is not zero, should input the CPO.

4. Panel Operation:

Turn on power on the back of the instrument and the instrument is standby. The following information will be shown on the panel.

S1	V6	T..... °C
%...	CP	000000

- 1) If spindle 1 is used at 6RPM, can directly press the start button on the panel to initiate testing. When display values become stable, read the CP value. (In order to make readout stable, if viscosity is relatively low, turn spindle about 3-5 rotations; if viscosity is relatively high, turn spindle about 1-2 rotations.)
- 2) If what the user chooses is not the default value S1 and V6, press spindle select and speed select button to enter selecting status, then press confirm button to confirm. When set-up is done, press start button to initiate testing.

NOTE: Spindle S0~S4 and V0.1/V0.3/V0.6/V1.5/V3/V6/V12/V30/V60 will scroll displaying repeatedly.

Explanations:

S: Stands for spindle number. S1 refers to spindle 1. (default value)

V: Stands for rotation speed. V6 refers to 6RPM. (default value)

T:..... °C: Temperature

%.....: Test value's percentage rate of the full value of the measuring range.

Cp: mPa·S Viscosity value

For example: If the viscosity of the liquid to be tested is estimated

3000mPa·S, it is recommended to choose the combination S2, V6 or S3,

V0.3.

For the correlations between the combination of spindle and rotation speed and the corresponding viscosity range, please refer to the below table.

NDJ-5S measuring range table

Rotor Rotating speed	0#	1#	2#	3#	4#
6 rotations	100	1000	5000	20000	100000
12 rotations	50	500	2500	10000	50000
30 rotations	20	200	1000	4000	20000
60 rotations	10	100	500	2000	10000


NDJ-9S measuring range table

Rotor Rotating speed	0#	1#	2#	3#	4#
0.1 rotations	6000	60000	300000	1200000	6000000
0.3 rotations	2000	20000	100000	400000	2000000
0.6 rotations	1000	10000	50000	200000	1000000
1.5 rotations	400	4000	20000	80000	400000
3 rotations	200	2000	10000	40000	200000
6 rotations	100	1000	5000	20000	100000
12 rotations	50	500	2500	10000	50000
30 rotations	20	200	1000	4000	20000
60 rotations	10	100	500	2000	10000

Attentions

1. When this viscometer works under room temperature, the tolerance of the test temperature should keep within $\pm 0.1^{\circ}\text{C}$. Otherwise measurement accuracy will be substantially affected.
2. Pay attention to test values and their percentage rates of The entire measurement range. When numbers are too high or too low, spindle or rotating speed should be changed to keep the percentage within 15%~85%. Otherwise measurement accuracy will be affected.
3. Viscometer should be used under allowed voltages and frequencies. Otherwise measurement accuracy will be affected.
4. Be careful when spindle is assembled or unassembled. Slightly lift up direction connector. Don't handle it with big force.
5. Keep spindle and direction connector clean after use.
6. When instrument goes down, it should be held by hands to avoid vibrations which could do damages to the pivot.
7. When instrument is moved or transported, direction connector should be covered by protection cap.
8. Many of suspension liquids, emulsions, high polymers and other high viscosity liquids are non-Newtonian liquids, whose viscosity values are changing with changes of shear rate and time. It is normal that the test results could be different on different spindles, rotating speed, or time.
9. In order to get accurate measurement, please pay attention to the following points:
 - 1). Control the temperature of the tested liquid accurately.
 - 2). Keep environmental temperature uniform.
 - 3). Keep spindle and liquid to be tested under constant temperature at the same time so as to keep the same temperature for both.
 - 4). Use bigger spindle with higher rotating speed for low viscosity liquid. Use smaller spindle with lower speed for high viscosity liquid.
 - 5). Keep the surface of the spindle clean.

No. 0 SPINDLE

Installation Procedure:

Screw on spindle 0 (5) ➡ put on fixing tube (2) ➡ screw on bolt (1) and get it fixed ➡ put on rotating sleeve 0 (4) ➡ screw bolt (3) ➡ tight and get it fixed.

No. 0 spindle table

Speed	60	30	12	6	3	1.5	0.6	0.3
Range	10	20	50	100	/			

Appendix: Spindle 0 installation procedure and usage.

- 1). Should consider instrument components operation when spindle 0 is used.
- 2). Pour 25~30ml liquid to be tested into lower sleeve.
- 3). Gently put outer testing tube from down to up on the fixing sleeve and get it fixed. (On the testing tube is a V-shaped groove)
- 4). The same method as the above applies when bottomless testing tube is used. Take V-shaped groove on the inner wall of outer testing tube as the immersed level line for spindle 0 when outer testing tube and spindle are immersed in the liquid.

MAINTENANCE

1. Clean the shell and internal components of the viscometer regularly, especially the measuring head, to avoid contamination and bacterial growth.
2. Do not use any acidic or alkaline cleaning agents.
3. Viscometers require frequent calibration to ensure accuracy and repeatability. It is recommended to calibrate before each use, while calibrating at different temperatures. (Please be calibrated by professional technicians.)
4. Pay attention to avoid collision and falling during use, especially the measuring head.

Packing List

Name	Quantity
Digital viscometer (Master device)	1
Lifting rod and slider	1
Protection Bracket	1
Spindles No.1,2,3,4	1 for each
Spindles No.0	1
Power cord	1
User Manual	1

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YH CONSULTING LIMITED.
C/O YH Consulting Limited Office 147, Centurion House, London Road, Staines-upon-Thames, Surrey, TW18 4AX

VEVOR[®]

TOUGH TOOLS, HALF PRICE

Technical Support and E-Warranty Certificate

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Importado a AUS: SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD NSW 2122 Australia
Importado a EE. UU.: Sanven Technology Ltd. Suite 250, 9166 Anaheim Place, Rancho Cucamonga, CA 91730
E-CrossStu GmbH