# **VEVOR<sup>®</sup>**

## **TOUGH TOOLS, HALF PRICE**

Technical Support and E-Warranty Certificate [www.vevor.com/support](http://www.vevor.com/support)

## **JADE STONE HEATING PAD**

**MODEL: H23Y1**

We continue to be committed to provide you tools with competitive price.

"Save Half", "Half Price" or any other similar expressions used by us only represents an estimate of savings you might benefit from buying certain tools with us compared to the major top brands and does not necessarily mean to cover all categories of tools offered by us. You are kindly reminded to verify carefully when you are placing an order with us if you are actually saving half in comparison with the top major brands.

# VEVOR®

TOUGH TOOLS, HALF PRICE

## JADE STONE HEATING PAD

MODEL: H23Y1

### NEED HELP? CONTACT US!

Have product questions? Need technical support? Please feel free to contact us:

**Technical Support and E-Warranty Certificate**  
[www.vevor.com/support](http://www.vevor.com/support)

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

	<p><b>Warning-</b>To reduce the risk of injury, user must read instructions manual carefully.</p>
	<p><b>CORRECT DISPOSAL</b></p> <p>This product is subject to the provision of European Directive 2012/19/EC. The symbol showing a wheellie bin crossed through indicates that the product requires separate refuse collection in the European Union. This applies to the product and all accessories marked with this symbol. Products marked as such may not be discarded with normal domestic waste, but must be taken to a collection point for recycling electrical and electronic devices</p>
	<p>Compliance is a EC &amp; UK security certification.</p>

**FCC INFORMATION**

**CAUTION:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment!

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This product may cause harmful interference.
- 2) This product must accept any interference received, including interference that may cause undesired operation.

**WARNING:** Changes or modifications to this product not expressly approved by the party responsible for compliance could void the user's authority to operate the product.

**Note:** This product has been tested and found to comply with the limits for a Class B digital device pursuant to Part 15 of the FCC Rules, These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This product generates, uses and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this product does cause harmful interference to radio or television reception, which can be determined by turning the product off and on, the user is encouraged to try to correct the interference by one or more of the following measures.

- Reorient or relocate the receiving antenna.
- Increase the distance between the product and receiver.
- Connect the product to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for assistance.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

## IMPORTANT SAFEGUARDS

Thank you for using this product. In order to make sure that you can operate the machine correctly, read this instruction carefully before operation and keep it properly for future reference. Please be sure to read the precautions and safety rules on this page to ensure your safe use.

This manual will outline safety warnings and precautions, operating, maintenance and cleaning. The warnings and instructions reviewed in this manual cannot cover all possible conditions and situations that may occur. Caution and common sense are not built into this product, since we believe that the users will comply with these codes. If you have any problems, please call us. Do not try to fix the problems yourself, or it will cause injury or property damage.

### **Please read ALL the instructions before using your machine.**

1. Please check the product before each use. If you find that the jade is broken or the product is damaged, please do not use it.
2. Do not let children, minors, and pets get close to or touch the machine.
3. Do not immerse the appliance in water or any other liquid.

4. Do not disassemble the product.
5. Using any accessories or attachments that are not provided or recommended by the manufacturer may cause injury.
6. Please start to use after the collect is fixed well.
7. Do not leave the product hanging over the edge of table or counter or let it come into contact with any hot surfaces.
8. Do not place the product on or near hot air, electric burners or heaters.
9. Do not apply this unit to any other purposes than the indicated usage.
10. KEEP IT OUT OF CHILDREN'S REACH.
11. DO NOT CLEAN IT WITH ANY ABRASIVE MATERIAL.
12. Please clean it regularly, at least once a month.
13. Do not cover the mat with a quilt or blanket or pillow and so on when using the mat.

**Note:** The quilt/blanket/pillow...may be burned when the heat of the mat is high.

14. Do not partially cover this pad to avoid partially overheating.
15. Do not use this pad on a latex/silicone/sponge bed mattress or sofa that is too soft.
16. Do not use this pad together with other heating products
17. Do not use while unconscious or while sleeping.
18. Do not use on an infant.
19. Do not use if you are pregnant or think you may be pregnant
20. Do not use over any malignancy. Anyone with diabetes, poor blood circulation or with disabilities should consult their physician before use.
21. Place pad directly where heat is desired.

**Note: Burns can occur regardless of control setting. Check skin under pad frequently to avoid burning and blistering.**

22. Do not use in any oxygen enriched environment or near equipment with stores or emit oxygen.
23. Do not fold pad and apply weight (Do not fold sharply).

24. Do not use pins,nails,screws,or any other metallic means to fasten this pad in place
25. Never pull the pad by the power supply cord or use the cord as a handle.
26. Carefully examine inner cover before each use, discard the pad if inner covering shows any signs of deterioration such as blistering or cracking.
27. If smoke is ever detected from this device, immediately remove the power cord from the outlet.
28. Do not use this pad with liniments,salves or ointments that contain heat producing ingredients, skin burns could result.
29. Use wherever hot applications are desirable for personal comfort, and whenever recommended by your physician for the relief or pain.
30. Keep you pad in its carrying bag or in a dry,cool place
31. Do not hang pad by the power cord.
32. Do not use this pad in the presence of flammable mixture.
- 33. LIFE SAFETY WARNING:** Do not use this product after drinking alcohol or taking other medications.
- 34. LIFE SAFETY WARNING:** Pregnant women, skin diseases, heart diseases, hypertension, brain diseases, and other diseases should consult a doctor before using this product.
35. This appliance is not intended for use by persons (including children) with educed physical, sensory or mental capabilities, or lack of experience or knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safe.
- CAUTION:**
- If you have any health concerns, do not hesitate to consult your doctor before using this product.
- Never leave the appliance unattended under any circumstance.Never use this product directly on swollen or inflamed areas of skin.
- Do not use this product if you suffer from a physical ailment that would limit or impair your ability to operate the controls. This unit should not be used by children without adult supervision.

## LIFE SAFETY WARNING

The list below contains known conditions when the far infrared heating pad should not be used or only be used under the supervision of a physician. This list is not exhaustive. If you suffer from any serious condition that might be affected by the use of the product, please consult with your physician before using the heating pad.

- Do not use directly over open wounds
- Do not use directly over acute injuries.
- Do not use directly over face or genitals.
- Do not use if you are insensitive to heat, including paraplegic and quadriplegic or have poor blood circulation that affects your ability to feel heat or cold sensations.
- Do not use if you are immobile or incapacitated or not able to remove the heating pad from your body.
- Do not use this product with infants or children under the age of 6. Children under age of 14 should always use the pad under adult supervision. In the case of pregnancy, diabetes, hemophilia, malignancy or other serious conditions, consult with your physician before applying the heating pad.

## BEFORE USING

Before using the appliance for the first time, peel off any promotional and packaging materials and check that the product is clean and dust-free. Please dispose of all plastic bags carefully and keep them away from children and pets. Check all components provided according to the list in this manual.

First, open the bag and confirm all components are included. Then perform the following:

Step1: Connect the heating pad to the digital controller.

Step2: Connect the Digital Controller to the AC-DC Adapter

Step3: Plug the Power Cord into the wall outlet

The Heating Pad is now assembled.

**Note: Before each use, it is necessary to check whether the product is intact and whether the connector links are reliable.**

## PRODUCT PARAMETERS

Model	H23Y1
Power adapter Input:	AC100-240V 50/60Hz
Power adapter Output:	DC24V 1.5A

**Note: Please use the correct power source voltage.**

## CONTROLLER INSTRUCTIONS

## **ON /OFF:**

The digital controller allows you to turn the unit “ON” or “OFF”, simply by pressing the ON/OFF button.

## **Vibration Massage:**

The digital controller has a vibration massage function, it has three vibration levels: level 1, level 2, level 3. By pressing the VIBRATE button, you can adjust the vibration from level 1 to level 3 or shut off the vibration. Each vibration will last for as long as 20 minutes at a time. Press the VIBRATE button for the first time, the vibration level is level 1. Press the VIBRATE button for the second time, the vibration level is level 2. Press the VIBRATE button for the third time, the vibration level is level 3. Press it again, the vibration is shutoff. To continue to use the vibration, Press the VIBRATE button again, it will work from level 1.

## **Heat:**

The digital controller allows you to adjust the heat to a temperature that you are comfortable with. Simply press the (TEMP+) button to increase the temperature or press the (TEMP-) button to lower the temperature by 1°F or 1°C increments, the temperature range is 103°F~159°F (40°C~70°C) .

When you press the TEMP+ or TEMP-, the TEMP. (°F or °C) light will be on, then starts blinking and the screen also starts blinking four times showing the temperature you set. Then the screen returns to show the current temperature of the pad. If the HEATING light is on all the time, it means the digital controller is heating up in order to reach the temperature that you set. If the HEATING light is blinking, it means the current temperature of the pad is close to the temperature that you set, the digital controller will slow down heating up. If the HEATING light is NOT on, it means the current temperature of the pad has already reached the temperature that you set so the digital controller has stopped heating up; Once the current temperature of the pad drops to the temperature settings, the digital controller will heat up and the HEATING light will be on again.

**Important Note:** This Far Infrared Heating Pad works a little differently than your standard coil heating pad. Please allow 5 to 10 minutes for the Far Infrared energy in the stones to build up and for the pad to heat up. Once you are ready to use, please adjust the temperature to your tolerance.

**Timer:**

The digital controller allows you to set the timer automatically to shut off. By pressing the (TIME+) or (TIME-) buttons you can adjust the time by 15 minutes increments. The timer can be set for as long as 240 minutes at a time. Once the timer counts itself all the way down, it will automatically shut off. When you press the TIME+ or TIME -, the TIME (Min) light will be on, then starts blinking and the screen also starts blinking four times showing the time you set. Then the screen returns to show the current temperature of the pad. If you do not adjust the timer, the unit will count up to 4hours and shut itself off.

**Important Note: Always unplug this device when not in use.**

## TIPS

1. After taking a shower, using this product will help the body relax more.
2. Before use, set the temperature to the lowest level and wait for the human body to adapt to the temperature before adjusting to the desired temperature.
3. Before each use, make sure that all connections are tight and snug.
4. Before each use, plug the pad into the wall outlet, the use of a surge protector is recommended.
5. Before each use, make sure that all of the functions of the digital controller are working correctly.
6. It is suggested that you get in a comfortable position, either sitting or lying down.
7. It is suggested that you place a small towel (folded in half) as a protective layer between you and the stones. There will not be a loss of effectiveness if the heat has to travel through the protective layer or even your clothing.
8. In other words, the stones do not need to be touching your skin to be effective. For safety reasons, the stones should not be touching your skin directly.
9. Depending on the thickness of the protective layers, you will need to adjust the temperature to your tolerance level. The temperature should never be set hotter than you can comfortably tolerate it

10. The pad should be positioned directly over the affected area. You may either lie on the pad, wrap the pad around the effected body part or lay the pad on your body or body part. Remember to use a protective layer (ex.Small towel) for both you and the pad to keep it clean.

11. When using this pad, please have the jade stones facing your body.

12. There is no specific protocol for any individual using a heating pad. With that being said everyone responds differently to different modalities and it is suggested that you limit the use of the pad to 20minutes per session of continuous use for the first several weeks to see how you tolerate the heat. Once you have become comfortable with the use of the pad, you can begin to use it for longer durations of continuous use. Please note that the longer the duration, the lower the heat needs to be, learning to adjust the heat to your tolerance will be a key to achieve optimal results.

## CLEANING

1. Before cleaning, the product plug must be unplugged. Put the product in a power-off state.

2. Use a soft and moist cloth for cleaning, do not rinse directly with water.

3. Please keep accessories and products together.

4. Do not fold excessively during storage to avoid damaging internal wiring.

- Do not submerge controller, Converter or pad in water.
- Do not dry clean any part of the heating pad.
- Do not use bleach or detergent as cleaning agents.
- Do not use a wringer
- Do not put in Microwave.
- Do not put in washing machine or dryer.
- Do not Iron
- Do not use mothballs or repellent sprays of any kind
- Do not put any of the power cords into water
- Do not use unless the pad is completely dry.
- Do not attempt to clean each individual jade stone.

## 5. Acceptable cleaning method

- Wipe the pad before the initial use and after each use.
- Use a gentle nonabrasive soap with a slightly moist towel.
- Let the heating pad dry completely before using.

## 6. Routine maintenance for safe performance

- Inspect the heating pad both before and after each use.
- Examine the surface of the heating pad for any sign of abnormal or excess wear and tear including blistering and cracking
- Examine the power cord for cracks, blisters, frays or exposed wires

## **STORAGE**

To store your heating pad, completely clean the unit, let it completely dry and place the heating pad back into the supplied storage bag. Store the heating pad in a cool and dry place away from any risk of damage while the pad is not in use

Only transport the heating Pad inside of the storage bag provided. Do not let the pad get wet or exposed to any moisture.

- Do not hang the heating pad by the power cord.
- Do not bend or put any additional weight on the heating pad.
- Do not allow the power cord to be pinched.

UTK TECHNOLOGY SHENZHEN LIMITED

**Address:** 501, A1 Buiding, Yinlong Industrial Park, No.292, Shenshan Road  
Longgang District Shenzhen, Guangdong CN 518118

**Manufacturer:** Shanghaimuxinmuyeyouxiangongsi

**Address:** Shuangchenglu 803nong11hao1602A-1609shi, baoshanqu, shanghai  
200000 CN.

**Imported to AUS:** SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD NSW 2122  
Australia

YH CONSULTING LIMITED.

C/O YH Consulting Limited Office 147,  
Centurion House, London Road,  
Staines-upon-Thames, Surrey, TW18 4AX

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Mainzer Landstr.69,  
60329 Frankfurt am Main.

**VEVOR<sup>®</sup>**

**TOUGH TOOLS, HALF PRICE**

**Technical Support and E-Warranty Certificate**

**[www.vevor.com/support](http://www.vevor.com/support)**

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**Importerad till AUS:** SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD NSW 2122
E-CrossStu GmbH