VEVOR

Affordable. Reliable. Home Improvement.


BODY SCULPTING MACHINE

MODEL:BCE1

VEVOR
Affordable. Reliable. Home Improvement.

Body Sculpting Machine

MODEL:BCE1

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

	Warning-To reduce the risk of injury, user must read instructions manual carefully.
	This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:(1)This device may not cause harmful interference, and (2)this device must accept any interference received, including interference that may cause undesired operation.
	This product is subject to the provision of European Directive 2012/19/EC. The symbol showing a wheelie bin crossed through indicates that the product requires separate refuse collection in the European Union. This applies to the product and all accessories marked with this symbol. Products marked as such may not be discarded with normal domestic waste, but must be taken to a collection point for recycling electrical and electronic devices

PRODUCT INTRODUCTION


This product falls under the category of personal care and massage devices. It generates suction through a

negative pressure pump, designed to provide massage and relieve mild muscle soreness.

The device comes with two detachable suction cups for flexible use.

Please note that the main unit is not waterproof.

It contains a built-in rechargeable lithium battery (3.7V / 1400mAh) and supports charging via a DC 5V/1A or higher power supply.

PRECAUTIONS

Please read the instructions carefully before using the

device.

1. Keep the device away from water.

2. Do not immerse or drop into the water or any other liquid.

3. If the device falls into water, unplug it immediately. Do not place your hand in the water.

4. A device should never be left unattended while it is plugged into the mains.

5. This device can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning usage of the appliance in a safe way and understand the hazards involved.

6. Use the device only for its intended purpose, as indicated in this manual. Do not use any accessory which is not recommended by the manufacturer.

7. Do not drop the device or insert any object in the openings.

8. Do not overuse, maximum 5 minutes in one area.

9. For people with loose muscles, this device may leave a bruised mark after usage. In case bruising occurs, stop usage immediately. Normally it will disappear gradually after 1 week.

10. If your skin is under nursing, or under any unhealthy condition, please consult your doctor before using the device.

11. Do not open the unit under any circumstances.

12. Never leave the unit unattended, especially when children are present

13. The device should only be used with the power supply supplied with the device.

OPERATION INSTRUCTION

Charging

1. When battery is not enough to last one cycle (10 minutes usage), all LED indications will flash 5 times at 0.5Hz frequency for warning, then, keep regular working mode for continuous usage.

2. During usage, when battery is too low to continue, the unit will power off after all LED indications flash 5 times at 0.5Hz frequency.

3. Plug the Type-C end of the cable into the charging port on the device.

4. Plug the USB end of the cable into a standard USB charging plug.

5. During charging, the LEDs on main unit will illuminate in sequence, and form a full circle after fully charged.

6. A full charge requires up to 2 ~ 4 hours and will provide up to 60 minutes of usage. Once fully charged, unplug the adaptor from the socket.


Instructions for usage


1. Choose the cup according to the size of the area to

be treated.

2. Gently place the suction cup to main unit at a 30 degree angle. Remember, the “pointing mark” on cup is to the lower right corner.

3. Twist the cup clockwise, until a “click” sound. “Pointing mark” is in the middle of the device.

4. Press Power button to start the device. Press once again to switch mode between Body Mode and Face mode.

5. Under Body Mode, press LIPO button to adjust suction level with related indication illuminated. The suction level will change in the following sequence by further pressing the LIPO button. At initial usage, start with level 1. Gradual- ly increase level as the body becomes accustomed to suction stimulation.

Level 1 → Level 2 → Level 3 → Suction OFF

6. Press EMS button to adjust EMS level with related indication illuminated. The EMS level will change in the following sequence by further pressing the EMS button. At initial usage, start with level 1. Gradually

increase level as the body becomes accustomed to EMS stimulation.

Level 1 → Level 2 → Level 3 → EMS OFF

7. Under Face mode, only EMS is activated, no suction generated. Use the device without either cup.

8. Once the settings are as desired, firmly place the cup against your body and keep a slow movement.

RECOMMENDATION: It is ADVISED to apply cream, lotion, body oil in targeted treatment area to have a smooth movement on the body and less electric shock feeling caused by EMS.

9. Treatment time is recommended as 2-3 minutes per area every day, maximum 30 minutes to cover the whole body.

CAUTION: Do NOT overuse, maximum 5 minutes in one area.

Note: It is normal for skin to turn red after usage, which will disappear after 1 day of non-usage.

10. During usage, keep a slow and continuous movement. Simply hold the device firmly against your body to prevent it from falling, no need to push it hardly.


CAUTION: Do NOT stay in one position for a long time.

11. When both LIPO and EMS are on, you should feel both suction of your skin and light electric shock which is normal phenomenon. If either becomes uncomfortable, press the RELEASE button for immediate stop. Adjust the level to a lower setting to

re-start the usage.

12. Once done, press Power button to turn off the device. The device includes an auto-shutoff after 10 minutes.

13. Simple twist cup anti-clockwise to remove from main unit.

CLEANING AND MAINTENANCE

Cleaning

1. Make sure the unit is un-plugged and is powered off.

2. Remove cup from the main unit.

3. Replace the old filter with a new one, as shown in the image.

Replacement frequency: once a month

Replacing filter when the suction power is less.

4. Wash clean the suction cup with tap water.

5. Lightly dampen a clean cloth to clean the main unit.

IMPORTANT: KEEP the unit away from water. Do not flash or immerse the unit with water.

IMPORTANT: Never use detergent to clean the main unit or suction, to avoid damage.

Storage

After usage, make sure to power off the unit and remove the adapter from the socket. Store the unit away from areas with high humidity or high temperature.

If it is not going to be used for a long time, thoroughly clean the main unit and the suction cup before storage.

Environment

This product should NOT be treated as household waste when you wish to dispose it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring the device is disposed correctly, you will help to reduce negative consequences to the environment

and human health, which may otherwise be caused by inappropriate disposal of the device. This recycling of materials will help to conserve natural resources.

BATTERY DISPOSAL

To minimize hazards to health and the environment at the end of this product's life，laws dealing with Waste Electrical and Electronic Equipment（WEEE）and The Waste Battery Directive require you to dispose of this product at a suitable collection facility where it will be sent in order to remove the batteries and for appropriate recycling.Please contact your local authorities for more details on recycling and safe disposal of these in your area.

TROUBLE SHOOTING

PROBLEM	PROBABLE CAUSE
Unable to turn on or off	The machine battery is dead.
Unable to charge	• The charger is not properly plugged in. • The charger being used does not match the correct specifications. • The charger or the data cable is faulty.
During use, the gear indicator light flashes, or after charging, it can only be used for a few minutes. The device suddenly stops working while in use.	• Low battery power • Automatic shutdown - The machine operates in a cycle of 10 minutes. After 10 minutes, the product will automatically shut down: • This may be due to a malfunction of the built-in rechargeable battery in the main unit.

Manufacturer: Shanghaimuxinmuyeyouxiangongsi

Address: Shuangchenglu 803nong11hao1602A-1609shi, baoshanqu,

shanghai 200000 CN.

Imported to AUS: SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD

NSW 2122 Australia

Imported to USA: Sanven Technology Ltd. Suite 250, 9166 Anaheim

Place, Rancho Cucamonga, CA 91730

UK REP YH CONSULTING LIMITED. C/O YH Consulting Limited Office 147, Centurion House, London Road, Staines-upon-Thames, Surrey, TW18 4AX

EC REP E-CrossStu GmbH
Mainzer Landstr.69,
60329 Frankfurt am Main.

VEVOR

Affordable. Reliable. Home Improvement.


MACHINE DE SCULPTURE CORPORELLE

MODÈLE : BCE1

VEVOR
Affordable. Reliable. Home Improvement.

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Importerad till Australien: SIHAO PTY LTD. 1 ROKEVA STREETEASTWOOD
Importerad till USA: Sanven Technology Ltd. Suite 250, 9166 Anaheim